



Software Tools and Services for chrono and psycho physiology

PRANA Viewer


Sleep Scoring Report

The PRANA® software for biosignal processing has convenient reporting capabilities.

The software report generator can configure and perform the instant edition and printing of comprehensive study reporting as computer files containing the relevant data and graphical representations.

The software continuous scoring report computes and displays the statistics related to the sleep-wake stages and contains:

- o the patient and study information, the hypnogram and general sleep parameters,
- o the internal sleep structure with stage characteristics, sleep-wake stages duration, latency, stability and cycle characteristics,
- o a set of graphics summarizing sleep-wake stages proportion, cumulated stage duration, etc.




Sleep scoring report generator



Sleep/wake stage scoring report

The software iterative scoring report for iterative vigilance tests (MSLT, MWT, others) computes and shows the temporal distribution of the sleep-wake stages (latency and proportion) as a graphic summary with the corresponding data available as text file.



| Region | Stage |
|--------|-------|--------|-------|--------|-------|--------|-------|--------|-------|--------|-------|
| 0 | Wake |
| 1 | Light |
| 2 | Deep |
| 3 | REM |

Iterative scoring report : graphical representation and corresponding data

Sleep Trend Report

The PRANA® software for biosignal processing has extensive reporting capabilities.

The software report generator can configure and perform the instant edition and printing of comprehensive study reporting as computer files containing the relevant data and graphical representations.

The software analyses report computes and displays the distribution by sleep stages and sleep cycles of any time series obtained after execution of a software optional plug-in. The report contains:

- o the time course of selected parameters and the concurrent hypnogram,
- o the mean value of selected parameters by sleep stage in sleep cycles and by sleep phase in sleep cycles,
- o the slope and intercept values of selected parameters by sleep phase in sleep cycles,
- o the cumulative value of selected parameters by sleep phase.



Sleep trend report generator



Sleep analysis report

Event Marking Report

The PRANA® software for biosignal processing has powerful reporting capabilities.

The software report generator can configure and perform the instant edition and printing of comprehensive study reporting as computer files containing the relevant data and graphical representations.

The software event marking report computes and shows the distribution of any type of event by sleep stages and sleep cycles and contains:

- o the patient and study information,
- o the global event type statistics (number and mean duration),
- o the internal sleep repartition and the corresponding statistics.



Event marking report